



# 2026 MD Glow

*Treats All Skincare Needs*

**Multifunctional Skin Care Platform 9 in 1**



FINE LINES + WRINKLES



ELASTICITY + FIRMNESS



EVEN TONE + VIBRANCY



SKIN TEXTURE



BROWN SPOTS



OILY + CONGESTED SKIN



ENLARGED PORES



**NINE DIFFERENT HANDLES**



Hydro Dermabrasion



Polymer Atomizing Spray



Ultrasound Import



Bipolar RF



Multi-pole Microcurrent



Bipolar Microcurrent



Cold Hammer



Ultrasonic Vibration




Mousse Bubbles

## DEEP CLEANING




### Hydro Dermabrasion


Open pores, soften the stratum corneum, deeply cleanse the skin, deeply moisturize, remove blackheads, improve the dark yellow skin, and make the skin clean, white and supple.



Solution is sprayed through the inlet of nozzle tip. It rotates clockwise.



The waste and remaining solution are sucked to container through inlets.



Solution are supplied to skin surface through nozzle tip.



Hard turbine hydro cleaning head



Microdermabrasion head



Soft silicone hydro cleaning head



Pipe cleaning head



## ***DEEP CLEANING***

### ***Mousse Bubbles Handle***

The handle produces dense foam containing a large amount of active oxygen components. Using the high oxidative properties of active oxygen, the hydroxide ions can quickly penetrate into the skin, so as to achieve the effect of cleaning, moisturizing and rejuvenating the skin.



### ***Ultrasonic Vibration Handle***

Remove dead skin by high-frequency vibration waves, remove blackheads, soften skin keratinocytes, and deeply cleanse skin dirt.



## ***MOISTURIZING***



### ***Ultrasound Import***

At a frequency of 1 million to 3 million vibrations per second to quickly open the pores, allowing the essence to penetrate into the basal layer of the skin.



### ***Polymer Atomizing Spray***

Combined with liquid products, the liquid is converted into mist-like small molecules by high-pressure atomization, and then quickly poured into the dermis by high-pressure spraying to achieve a deep hydrating effect.



## ***TIGHTENING & LIFTING***



### **Microcurrent Handles**

Microcurrent uses low-voltage electricity to stimulate muscle, adenosine triphosphate (ATP) cell growth, and collagen development in the dermis on the face.



### **Biopolar RF Handle**

RF treatment is a non-invasive way to achieve skin tightening and make you look younger. The RF energy will heat the deep skin layers (dermis) of your body to produce collagen.





***Cold & Warm Hammer Handle***

Cold Hammer 0°C: Ice skin, shrink pores, calm skin

Warm Hammer 45°C: increase the temperature of the skin surface, accelerate blood circulation, and increase the absorption rate when used with products.



***SCOPE OF TREATMENT***

**Deep Cleaning**

Remove Blackheads, Mites Clear, Remove grease dust,  
Remove Dead Skin


**Skin Rejuvenation**

Skin Moisturizing, Shrink Pores, Essence Import, Skin  
Whitening, Acne Treatment

**Skin Tightening**

Face Lifting, Eye Lifting, Wrinkle removal





## ***SPECIFICATION***

Power	600w
RF frequency	3Mhz
Screen	8.4inch
Vacuum	0-90kp
Flow Capacity	30L Max
Probe of dermabrasion	11pics
Voltage	100V-240V 50Hz-60Hz



7696 183A Toll Road, Suite 1C  
Leander, TX. 78641  
888-446-5727  
[www.MDLasers.com](http://www.MDLasers.com)